


i400 – Clinical excellence meets design-led style



The i400 combines the most advanced iPulse™ technology with an exclusive design for premium skincare establishments. This latest addition to the celebrated iPulse™ range perfectly combines aesthetics and clinical excellence.

i400



Treat large areas quicker

Like all iPulse™ systems, the i400 is designed to produce the most effective clinical results in the quickest time, at the lowest possible cost to the practitioner. The system's large 8.9 cm² beam of light means areas of the body are treated quicker with the i400, and the perfect level of energy generated by icePulse, a proprietary power supply, makes each pulse of light more effective, resulting in fewer treatments. This saves the practitioner time and money.

Lighter and sleeker

Cosmetic devices using light energy have traditionally been heavy, bulky machines, while smaller devices (sometimes known as mini-IPLs) compromise on clinical results. The advanced icePulse technology behind the i400 allows it to deliver the highest clinical performance from a small 31.8 lb. unit. This compact unit is portable within the treatment environment and fits neatly onto a mobile cart or standard desktop.

Two options

The i400 is available in two versions. One is a small table-top unit and the other comes with an integral cart designed to suit a clinical environment where the device needs to be used in different areas.



i400 – an IPL for the style conscious




Technical Specification

| | |
|---------------------------------|--|
| Light source: | twin Xenon lamps |
| Wavelength range: | 530 – 1200 nm |
| Spot size on tissue: | 8.9 cm ² |
| Pulse duration: | 10 – 110 ms |
| Pulse sequence: | single or triple pulse flash |
| Triple pulse: | 5 – 20 ms on-time 7 – 30 ms off-time |
| Repetition rate: | average of 0.25 Hz |
| Energy density (fluence) range: | 0 – 20 J/cm ² |
| Cooling: | fan, air-cooling |
| Dimensions for table-top: | 7.1" (high) 19.7" (wide) 15.6" (deep) |
| Dimensions for cart version: | 37.8" (high) 19.7" (wide) 15.6" (deep) |
| Weight: | 31.8 lbs. (table-top) 43.0 lbs. (with cart) |

Improved User-Interface

Adjustment of settings has been simplified so the device can be easily set up for different clients. The practitioner selects the treatment parameters based on a combination of skin types and desired results to provide exactly the right treatment for each individual client.



Touch sensitive control and LCD display allow for easy set-up.

Ready to use

The i400 plugs into a standard wall outlet and is ready to use immediately. Lamp cartridges provide the optimum light energy for treatment throughout their fixed lifetime; they are just as powerful on their 10,000th shot as their first. The computer alerts the operator when the lamp needs to be replaced...a ten second task for the operator just like changing a normal light bulb.



The i400 is available in two versions—table-top and standalone with integral cart.

All Inclusive Package

The i400 system features:

- One multi-procedural hand piece
- One lamp assembly
- Two pairs of safety glasses
- One user's manual
- One year parts and labor warranty
- FDA cleared for multiple indications

